

IoT Based Heart Monitoring System Using ECG

Karingala Suraj Raj - 180102026

November 2020


1 Main Objective

In the present world where people have become so busy with their work and with increased pollution, irregular eating habits most of the people are neglecting their health and have no time for health checkups. This could lead to serious health issues. So, the main objective of this project is to continuously monitor the heart of the people especially old people who either live in remote areas or have no one to look after them. The traditional 12 electrode ECG test is expensive and not portable, this project aims at providing cheap portable and continuous monitoring of heart without compromising much on the results. The data collected via sensors is sent to the cloud so that the healthcare personnel can access and analyse the data and can call the person for a checkup if they find any abnormalities.

2 Implemented Attributes


- AD8232 ECG Sensor is used to get the ecg data from the patient via 3 electrodes and to use it for analysis.
- Used the Blynk mobile application to display the ecg signal and BPM of the patient.
- Used Ubidots, a IOT platform as a cloud server to collect data from nodemcu over internet and display the information on its dashboard and also provide data for further analysis through MQTT protocol.
- An alert email and sms is sent to the patient and doctor from the ubidots cloud whenever there is any abnormality found or to remind the doctor's appointment or to remind the patient to take medicine on time.
- An LCD display displays the patients heart status regularly at server end.

3 Configuration Diagram



4 Sample Outputs

- The real-time raw data BPM will be calculated continuously and displayed.
- On pressing the Detailed analysis signal, few samples will get collected analysed and resulting BPM would be displayed in the app.



Sample output of the ECG Signal plotted on Serial Plotter.

5 Codes

The codes are attached with the document.

Code for Client (NodeMCU) :

```
1 //CODE FOR CLIENT NODEMCU
2
3 #include <ESP8266WiFi.h>
4 #include <PubSubClient.h>
5 #define BLYNK_PRINT Serial    // Comment this out to disable prints and
    save space
6 #include <SPI.h>
7 #include <BlynkSimpleEsp8266.h>
8
9 #define WIFISSID "Karingala"    // Put your
    WifissID here
10 #define PASSWORD "karingala@3579"    // Put
    your wifi password here
11 #define TOKEN "BBFF-kfFt6GXIEVhTrDQeK3yqdzxJJUSq0X"    // Put your
    Ubidots' TOKEN
12 #define MQTT_CLIENT_NAME "myecgsensor"    // MQTT
    client Name, please enter your own 8-12 alphanumeric character ASCII
    string;
13
    //it should be
    a random and unique ascii string and different from all other devices
14 char auth[] = "I2q95uptbB_wj5dXURbpagVw6yWmRSzd"; //Enter the Auth code
    which was send by Blink
15
16 #define VARIABLE_LABEL "ECG_graph" // Assigning the variable label, this
    will automatically create a variable in UBIDOTS Website
17 #define VARIABLE_LABEL_1 "BPM_graph" // Assigning the variable label, this
    will automatically create a variable in UBIDOTS Website
18 #define DEVICE_LABEL "ECG_Monitor" // Assigning the device label as
    created in UBIDOTS Website
19
20
21 #define SENSOR A0 // Set the A0 as SENSOR
22
23 char mqttBroker[] = "industrial.api.ubidots.com";
24 char payload[100];
25 char topic[150];    //payload and topic are for ECG_graph payload1 and
    topic1 are for BPM_graph
26 char payload1[100];
27 char topic1[150];
28 // Space to store values to send
29 char str_sensor[10];
30 char str_sensor1[10];
31
32 int bpm=0;
33 unsigned long oldMillis=0;
34 unsigned long newMillis=0;
35 long total=0;
36 long count=0;
```

```

37 long BPM=0;
38 float myecg=0.0;
39
40 WiFiClient ubidots;
41 PubSubClient client(ubidots);
42
43 void callback(char* topic, byte* payload, unsigned int length) {
44     char p[length + 1];
45     memcpy(p, payload, length);
46     p[length] = NULL;
47     Serial.write(payload, length);
48     Serial.println(topic);
49 }
50
51 void reconnect() {
52     // Loop until we're reconnected
53     while (!client.connected()) {
54         Serial.println("Attempting MQTT connection...");
55
56         // Attempt to connect
57         if (client.connect(MQTT_CLIENT_NAME, TOKEN, "")) {
58             Serial.println("Connected");
59         } else {
60             Serial.print("Failed, rc=");
61             Serial.print(client.state());
62             Serial.println(" try again in 2 seconds");
63             // Wait 2 seconds before retrying
64             delay(2000);
65         }
66     }
67 }
68
69 BlynkTimer timer;
70 void sendSensor()
71 {
72
73
74     // You can send any value at any time.
75     // Please don't send more that 10 values per second.
76     Blynk.virtualWrite(V5, myecg); //V5 is a virtual pin for ecg assigned
    in BLYNK app
77     Blynk.virtualWrite(V0, BPM);    //V0 is a virtual pin for BPM assigned
    in BLYNK app
78 }
79
80
81 void setup() {
82     Serial.begin(115200);
83     WiFi.begin(WIFISSID, PASSWORD);
84     // Assign the pin as INPUT
85     pinMode(SENSOR, INPUT);
86
87     Serial.println();
88     Serial.print("Waiting for WiFi...");

```

```

89
90 while (WiFi.status() != WL_CONNECTED) {
91     Serial.print(".");
92     delay(500);
93 }
94
95 Serial.println("");
96 Serial.println("WiFi Connected");
97 Serial.println("IP address: ");
98 Serial.println(WiFi.localIP());
99 client.setServer(mqttBroker, 1883);
100 client.setCallback(callback);
101 newMillis=millis();
102 Blynk.begin(auth, WIFISSID, PASSWORD);
103
104
105
106 // Setup a function to be called every 10 milliseconds
107 timer.setInterval(10L, sendSensor);
108
109 }
110
111 void loop() {
112     if (!client.connected()) {
113         reconnect();
114     }
115
116     sprintf(topic, "%s%s", "/v1.6/devices/", DEVICE_LABEL);
117     sprintf(topic1, "%s%s", "/v1.6/devices/", DEVICE_LABEL);
118     sprintf(payload, "%s", ""); // Cleans the payload
119     sprintf(payload, "{\\\"%s\\\":", VARIABLE_LABEL); // Adds the variable label
120     sprintf(payload1, "%s", ""); // Cleans the payload
121     sprintf(payload1, "{\\\"%s\\\":", VARIABLE_LABEL_1); // Adds the variable
        label
122
123
124 myecg = analogRead(SENSOR);
125 Serial.println(myecg);
126
127 //code to calculate BPM
128 if(myecg>500.0){
129     oldMillis=newMillis;
130     newMillis=millis();
131     long diff=(newMillis-oldMillis);
132     bpm=60000/diff;
133     total+=bpm;
134 }
135
136 count++;
137 if(count==500){
138     BPM=total/500;
139     Serial.println(BPM);
140     count=0;
141     total=0;

```

```

142 }
143
144 /* 4 is minimum width, 2 is precision; float value is copied onto
    str_sensor*/
145 dtostrf(myecg, 4, 2, str_sensor);
146 dtostrf(BPM, 2, 1, str_sensor1);
147 sprintf(payload, "%s {\"value\": %s}", payload, str_sensor); // Adds
    the value
148 sprintf(payload1, "%s {\"value\": %s}", payload1, str_sensor1); // Adds
    the value
149 Serial.println("Publishing data to Ubidots Cloud");
150 client.publish(topic, payload);
151 client.publish(topic1, payload1);
152 client.loop();
153 delay(10);
154 Blynk.run(); // Initiates Blynk
155 timer.run(); // Initiates SimpleTimer
156 }

```

Code for Server (NodeMCU) :

```

1 //CODE FOR SERVER SIDE NODEMCU
2
3 #include "UbidotsESPMQTT.h"
4 #define BLYNK_PRINT Serial
5 #include <Wire.h>
6 #include <LiquidCrystal_I2C.h>
7 #include <LinkedList.h>
8 // Set the LCD address to 0x27 for a 16 chars and 2 line display
9 LiquidCrystal_I2C lcd(0x27, 16, 2);
10
11 #include <ESP8266WiFi.h>
12 #include <BlynkSimpleEsp8266.h>
13 char auth[] = "I2q95uptbB_wj5dXURbpagVw6yWmRSzd"; //Enter the Auth code
    which was send by Blynk to your registered email id
14 int bval=0;
15
16 #define TOKEN "BBFF-i0MP7CIRKgFSC1vRMnp7g2cyK5DHsF" // Your Ubidots TOKEN
17 #define WIFINAME "Karingala" //Your SSID
18 #define WIFIPASS "karingala@3579" // Your Wifi Pass
19 #define VARIABLE_LABEL "ecg_graph" // Assinging the variable label
20 #define DEVICE_LABEL "ECG_Monitor" // Assinging the device label
21
22 Ubidots client(TOKEN);
23
24
25 int bpm=0;
26 unsigned long oldMillis=0;
27 unsigned long newMillis=0;
28 unsigned long startMillis=0;
29 long total=0;
30 long BPM=0;
31 String s="";
32 int a;

```

```

33
34 void callback(char* topic, byte* payload, unsigned int length) {
35     Serial.print("Message arrived ");
36     Serial.print(topic);
37     Serial.print("] ");
38     for (int i=0;i<length;i++) {
39         s+=(char)payload[i];
40         Serial.print((char)payload[i]);
41
42
43     }
44     Serial.println();
45     a=s.toInt();
46     s="";
47     //Serial.println(a);
48 }
49
50 BLYNK_WRITE(V1){
51     int pinvalue = param.asInt();
52     Serial.println(pinvalue); //pinvalue is the instruction that we are
53                               //getting from a button in BLYNK app
54     bval=pinvalue;
55 }
56 BlynkTimer timer;
57 void sendSensor()
58 {
59
60
61     // You can send any value at any time.
62     // Please don't send more that 10 values per second.
63     Blynk.virtualWrite(V2, BPM); //V2 is a virtual pin assigned in BLYNK app
64     //to show BPM after detailed analysis in a gauge variable
65 }
66
67 void setup() {
68     // initialize the LCD
69     lcd.begin();
70
71     // Turn on the backlight and print a message.
72     lcd.backlight();
73     lcd.clear();
74     // Printing welcome message
75     lcd.print("Bonjoure!!");
76
77
78     // put your setup code here, to run once:
79     client.ubidotsSetBroker("industrial.api.ubidots.com"); // Sets the
80     //broker properly for the Industrial account
81     client.setDebug(false); // Pass a true or false bool value to activate
82     //debug messages
83     Serial.begin(115200);
84     client.wifiConnection(WIFINAME, WIFIPASS);

```

```

83  client.begin(callback);
84  client.ubidotsSubscribe(DEVICE_LABEL, VARIABLE_LABEL); //Insert the
    dataSource and Variable's Labels
85  Blynk.begin(auth, WIFINAME, WIFIPASS);
86  lcd.clear();
87  newMillis=millis();
88  startMillis=millis();
89  // Setup a function to be called every 2 seconds
90  timer.setInterval(2000L, sendSensor);
91  }
92
93  int count=0;
94  int check=0;
95  int result=0;
96  int flag=0;
97  int i=0;
98  float mylist[200]; // A list to store ecg values that could be sent for
    further detailed analysis
99
100
101 unsigned long DateMillis=0; //give the date of the arrival of doctor here
    in millis
102 unsigned long TimeMillis=0; //give the time when patient has to take
    medicine here in millis
103 int doc=0;
104 int med=0;
105
106 //code for doctor arrival
107 void doctorarrival(){
108     unsigned long cur=millis();
109     if(DateMillis-cur==0){
110         doc=1;
111         return;
112     }
113 }
114
115 //code for medicine time
116 void medicinetime(){
117     unsigned long cur=millis();
118     if(TimeMillis-cur==0){
119         med=1;
120         return;
121     }
122 }
123 void loop() {
124     // put your main code here, to run repeatedly:
125     doctorarrival();
126     medicinetime();
127     if(!client.connected()){
128         client.reconnect();
129         if(bval==1){
130             client.ubidotsSubscribe(DEVICE_LABEL, VARIABLE_LABEL); //Insert the
                dataSource and Variable's Labels
131             mylist[i]=a;

```



```

132     i++;
133     if(i>=200){
134         bval=0;
135         Serial.println(bval);
136         check=1; //A flag to tell the compiler to send the list for
detailed analysis
137     }
138     }
139     }
140     client.loop();
141     Blynk.run(); // Initiates Blynk
142     //code to display the patient's on lcd display when the results are
being computed.
143     if(bval==1 || check==1){
144         lcd.setCursor(0,0);
145         lcd.print("Monitoring.....");
146         delay(4000);
147         result=1;
148         bval=0;
149         check=0;
150     }
151     // For simplicity i am applyong the same formula as used for calculating
real time BPM for raw data on the list of collected ecg values
152     //we can send this list to MATLAB for further processing and come back
to nodemcu with the processed results.
153     if(check==1){
154         for(int j=0;i<200;i++){
155             if(mylist[i]>500.0){
156                 oldMillis=newMillis;
157                 newMillis=millis();
158                 long diff=(newMillis-oldMillis);
159                 bpm=60000/diff;
160                 total+=bpm;
161                 count++;
162             }
163         }
164     }
165     BPM=total/count;
166     Serial.println(BPM);
167     count=0;
168     total=0;
169     check=0;
170     result=1;
171 }
172
173 //code to display the final results of patient's heart condition on lcd
display
174 if(result==1 && flag==0){
175     lcd.clear();
176     lcd.print("BPM : ");
177     lcd.print(BPM);
178     //There could be many conditions applied afrom analysis but for
simplicity im just considering BPM less than and greater then 100
179     if(BPM<=100){

```

```

180     lcd.setCursor(0,1);
181     lcd.print("Status : OK");
182 }
183 else{
184     lcd.setCursor(0,1);
185     lcd.print("Status : CONSULT");
186 }
187 flag++;
188 }
189 timer.run();
190 }

```

6 User Manual

- Connect nodemcu to the LCD display at the server side and Connect nodemcu to the AD8232 sensor at the client side.
- Place the electrodes in a triangular formation and make sure you wear some footwear while performing the experiment.
- Make sure that you are not charging your laptop while performing the experiment as it could cause 50-60Hz AC noise. If possible try to power the client side setup using a battery.
- Adjust the brightness of the LCD display by rotating the potentiometer attached to the I2C chip with a screwdriver.
- Create an account on Ubidots Website and you will be mailed the Ubidots token to your email id which you need to give in the code.
- Install Blynk app from playstore and add gauges,buttons and lineplots as mentioned in the pic([Click Here](#)).
- Finally upload the corresponding codes to the corresponding Nodemcus and install the required libraries if they are not installed on your PC.
- To get sms and email notifications add events in the data section of ubidots website.

7 Video Demo and other links

- Drive link to the video demo - [Click Here](#)
- Code for Server Nodemcu - [Click Here](#)
- Code for Client Nodemcu - [Click Here](#)